

Sida loo qaato suufka dareeraha afka

1

2 x suuf buluug ah oo ku jira tuubo cad oo balaastiig ah

2 x bac cad oo balaastiig ah, oo ku calaamadeysan "VRD" iyo "PHL"

Kartoon ku calaamadaysan UKHSA URGENT MEASLES TESTING

Waxyaabaha Ku Jira

Foomka codsiga

Baqshadda cad ee balaastiigga ah

2

2 mins

X2

Baaketka ka soo saar suufka buluugga ah. Suufka buluugga ah ku xoq xariiqda cirridka ilaa 2 daqiiqo.

3

4

X2

Ku celi tallaabooyinka 2 iyo 3

Suufka qoyan ku celi gudaha tuubada cad oo bedel furka. Ku qor magaca, taariikhda dhalashada (DOB), iyo taariikhda maanta bukaanka tuubo kasta.

5

VRD

Tuubada ugu horeysa ku rid gudaha bacda cad ee VRD. K saar hawada, ka fuji xariijimaha cas ee kor kaga dhegan kadibna xir.

6

PHL

VRD

Tuubada labaad ku rid bacda cad ee ku calaamadaysan PHL. K saar hawada, ka fuji xariijimaha cas ee kor kaga dhegan kadibna xir.

7

UKHSA URGENT MEASLES TEST

PHL

VRD

Calaamadeey oo buuxi foomka codsiga oo labada bacood ee cad oo ay ku jiraan muunadaha waxaad ku ridaa kartoonka kadibna xir.

8

UKHSA URGENT MEASLES TESTING

UKHSA URGENT MEASLES TESTING

Kartoonka ku rid baqshadda cad, xir oo ku soo celi xirfadlaha daryeelka caafimaadka ama u soo dhiib baqshad qaadaha.

© Crown copyright 2023. Guide to oral fluid swabs for urgent measles testing. UMTOF: SOM Somali. UKHSA GW2022-838. Hagahan ayaa lagu soo dejiisan karaa luuqadaha kala ah Carabi, Bengaali, Boortaqis Braazil, Faransiis, Hindi, Pashto, Polish, Romaniyan, Ruush, Soomaali, isbaanish, Turki, Yukreeniyaan, Urduu, Yiddish, Yoruba.